

PLASTERLESS ARTICULATOR

Study Model Articulator

“Clamp-on” plasterless articulator eliminates mess, speeds mounting of any pre-cast study model!

Avoid the time-consuming, messy and material intensive task of mixing plaster for use with traditional articulators. Our Plasterless Articulator makes it fast and easy to mount dental casts—with no plaster needed!

Use it to expedite:

- Diagnosis.
- Case presentations.
- Splint fabrications.
- Ortho treatment and equilibration planning.
- Orthognathic surgery prep.

Securely mount and remove any size casts with the twist of a thumb-screw. Sliding, self-centering arm provides full anterior-posterior positioning and 360° rotation of upper cast. Quick-release hinge allows separation of mounted models. Solid steel construction with white enamel finish.

One complete articulator per box.

11510 Plasterless Articulator

