

# TACK-EASE™

## Disposable Tacking Tip

Tack-Ease™ instantly converts most curing lights into a tacking light, saves time, ends tip build-up.

*Convert your curing light into a tacking light*

Simply place a Tack-Ease Disposable Tacking Tip onto your existing curing light tip or lens and you have a 1.5mm diameter tacking light tip perfect for tacking veneers and indirect inlays or onlays.

Invented by a dentist, Tack-Ease prevents marginal excess luting cement from curing to help speed clean-up time. It prevents composite build-up on curing tips too! Semi-flexible design helps you apply gentle, controlled pressure and ensures proper seating. Light-blocking PVC construction provides a snug, non-slip fit on most tips and lenses.

50 disposable tips per box.

**11410** Tack-Ease Disposable Tacking Tips. Specify 8mm, 9mm, 11mm or 13mm tip size when ordering.

